

Vacuum ultraviolet (VUV) spectroscopy: SOXMOS

C.Suzuki, H. Funaba, N. Tamura

e-mail: csuzuki@nifs.ac.jp

1. Purpose / Application

Measurement of temporal evolution of vacuum ultraviolet (VUV) and soft X-ray (SXR) emission spectrum and impurity line intensities

2. Name of analysis (Kaiseki) data / module of MyView2

'soxmos': Uncalibrated (rough) wavelength vs. counts data for channels 1 and 2

'soxmos1': Calibrated wavelength vs. counts data for channel 1

'soxmos2': Calibrated wavelength vs. counts data for channel 2

3. General Description (Port, field line, time resolution, spatial resolution, number of channels, etc.)

3.1 Features

- 2 m Rowland circle mounting
- 88 degree grazing incidence
- Dual multichannel plate (MCP) detectors
- 600 g/mm or 133.6 g/mm interchangeable gratings
- Overall wavelength coverage : 2–130 nm
- Spectral resolution : ~0.01 nm (600 g/mm grating) or ~0.06 nm (133.6 g/mm grating)
- Minimum readout time for a full spectral portion : 16 ms

3.2 Apparatus

- An entrance slit, a grating, and two microchannel plates are arranged along a Rowland circle (see Fig. 1).
- Spectra on the microchannel plates are converted to visible images via phosphor screen, and then transmitted and bundled via a fiber optic conduit to a 2048 channel photodiode array detector (see Fig. 2).
- Two different spectral regions are simultaneously measured.
- Groove density of the grating is selectable from 600 and 133.6 g/mm.




Fig. 1. Inside view of the SOXMOS spectrometer body.




Fig. 2. Schematic view of the SOXMOS spectrometer system.

3.3 Arrangement

- The device has been installed at 7O port (AL01-06).
- VUV spectra are measured along a line of sight in a horizontally elongated cross section.
- Angle of the line of sight can be tilted from -2 degree to 2 degree.




Fig. 3. Arrangement of the SOXMOS spectrometer.

4. Requirement in use

The required wavelength range should be notified in advance.

5. Description of analysis (Kaiseki) data / module of MyView2

'soxmos': Uncalibrated (rough) wavelength vs. counts data for channels 1 and 2

'soxmos1': Calibrated wavelength vs. counts data for channel 1

'soxmos2': Calibrated wavelength vs. counts data for channel 2

6. Others

Wavelength calibration is required every time the wavelength range changes.

References

[1] J. L. Schwob, A. W. Wouters, S. Suckewer, and M. Finkenthal, Rev. Sci. Instrum. 58 (1987) 1601.